

UNIVERSITY OF WATERLOO

Physics 360/460 – Experiment 4 **THERMIONIC EMISSION OF ELECTRONS**

Purpose: To verify the form of the Richardson-Dushman equation of thermionic emission using a GRD7 tube and obtain a value for the thermionic work function for the tungsten filament.

References:

The following references (electronic copies provided on LEARN) detail the background theory for thermionic emission.

1. T. B. Brown, *Electronics*, New York: Wiley 1954, pg.216-222.
2. J. Millman, *Vacuum tube and Semi-conductor electronics*, New York: McGraw-Hill 1958 pg. 60-79.
3. R. L. Sproull, *Modern Physics*, 2nd ed., New York: Wiley 1956, pg. 434-443.


Procedure:

1. With anode voltage, $V_a = 300$ V, set the filament current, I_f , to give an anode current, $I_a(V_a = 300$ V) = 8.00 mA. Leave the filament current I_f constant and record the magnitude of the filament voltage, V_f .
2. Measure I_a vs. V_a for $V_a = 300$ V, 250 V, 200 V, 150 V and 100 V.
3. Repeat steps 1 and 2 for $I_a(V_a = 300$ V) = 4.00, 2.00, 0.8 and 0.4 mA


Note: The effective anode voltage is that between the anode and the midpoint of the filament. Thus it is necessary to subtract half of V_f (see circuit) to correct each value of V_a . Note the V_f will be different for each initial setting for I_a .

Apparatus:

Circuit Diagram



Electrode Structure:



Suggested Instrumentation:

- Ferranti GRD7 Guard Ring Diode
- Healthkit Power Supply Model IP-17 (Plate voltage)
- Taylor Multimeter (Plate current)
- Simpson 260 meter (Filament voltage)
- HP 3468A Multimeter (Filament current)
- Harrison Power Supply Model 6281A (Filament current)

Note: **Maximum Filament voltage: 6.5 volts**
Filament Dimensions: Effective length 14.5 mm x 0.130 mm diameter

Analysis: Plot the I_a vs. V_a to find the values of I_o for different filament currents I_f . From the values of I_o and T , make a plot of I_o/T^2 versus $1/T$ on a semi-log scale. Using this plot, determine the work function of the filament, w_o , in eV. Compare the your value of w_o to literature.

The following empirical equation can be used to relate the filament current to the temperature of the filament in degrees Kelvin

$$T = 60.2 \sqrt{B(1 + 83 * 10^{-6} B)} \text{ K} \quad \text{where} \quad B = \frac{I_f}{d^2} \frac{\text{A}}{\text{cm}^2}$$

Additional background information:

Thermionic Emission: Determination of the Work function

Richardson - Dushman Equation:

$$J_o = AT^2 e^{-w_o/KT} \quad (1)$$

J_o – current density

T – temperature of filament

w_o – is the work function at zero E - field

K – Boltzman constant

A – a constant

The Schottky effect accounts for the fact that an electric field acting on the surface of the metal will reduce the work function. The new work function is now:

$$w = w_o - \sqrt{e^3 E / 4\pi\epsilon_o} \quad (2)$$

w – work function at some value of the Electric field

E – Applied Electric field

e – electronic charge

Now the R-D equation, (1), becomes:

$$J = AT^2 e^{-w/KT}, \quad (3)$$

which can be alternately written as,

$$J = AT^2 e^{-w_o/KT} e^{\frac{\sqrt{e^3 E / 4\pi\epsilon_o}}{KT}}, \quad (4)$$

or

$$J = J_0 e^{\frac{\sqrt{e^3 E / 4\pi\epsilon_0}}{KT}} \quad (5).$$

Replacing the current density, J , with the current, I , and taking the natural log of both sides of (5) gives,

$$\ln(I) = \ln(I_0) + \frac{\sqrt{e^3 / 4\pi\epsilon_0}}{KT} \sqrt{E}, \quad (6)$$

which in terms of the anode voltage, V , is,

$$\ln(I) = \ln(I_0) + C\sqrt{V}, \quad (7)$$

From, (7) it is apparent that a plot of $\ln(I)$ vs. \sqrt{V} will allow a determination I_0 .

Finally, eqn. (1) (with J_0 replaced by I_0) can be used to determine the work function from a plot of semi-log plot of (I_0/T^2) vs. $1/T$.